OneTouch® Customer Care Line/ Ligne InfoSoins OneTouch®

1 800 663-5521

www.OneTouch.ca

The system described herein is covered by one or more of the following U.S. patents: 5,708,247, 5,951,836, 6,045,567, 6,156,051, 6,197,040, 6,241,862, 6,284,125, and D428,150. Use of the monitoring device included herein is protected under one or more of the following U.S. patents: 6,413,410, 6,733,655. Purchase of this device does not act to grant a use license under these patents. Such a license is granted only when the device is used with OneTouch® Ultra® Test Strips. No test strip supplier other than LifeScan is authorized to grant such a license. The accuracy of results generated with LifeScan meters using test strips manufactured by anyone other than LifeScan has not been evaluated by LifeScan.

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Additional Instructions:

Getting a Blood Sample from an Alternate Site

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Forearm and palm sampling means more options

Meters that use OneTouch® Ultra® Test Strips allow you to sample blood from your forearm or palm as well as your fingertip. Having these options can mean having more flexibility—and more comfort.

Sampling from your palm or forearm allows you to use your fingertips less often. These alternate areas have fewer nerve endings than your fingertip.

Getting a blood sample from your forearm or palm is different than the method for getting a sample from your fingertips described in your Owner's Booklet. At times, results obtained at the forearm or palm may vary from a fingertip measurement. Be sure to read the Owner's Booklet as well as these instructions before testing.

TABLE OF CONTENTS

Important information about sampling from your forearm or palm	2
Choose the right sampling site at the right time	2

Sampling instructions	4
Choosing your sample site	5
Preparing the OneTouch® UltraSoft® Blood Sampler	6
Getting a blood sample	10
Completing your test	12

Important information about sampling from your forearm or palm

Choose the right sampling site at the right time

Glucose levels in forearm or palm blood samples may differ significantly from levels in fingertip samples, especially when your blood glucose level is changing rapidly. During these times, forearm or palm samples may show these changes more slowly. Talk to your healthcare professional before you begin using your forearm or palm for sampling.

If you are testing:	Use blood sample from your:
Routinely before meals	Fingertip, forearm, or palm
Prior to or more than two hours after: <ul data-bbox="62 693 455 831" style="list-style-type: none">• a meal• a rapid-acting insulin injection or insulin pump bolus• exercise	Fingertip, forearm, or palm
When your blood glucose is changing rapidly, such as: <ul data-bbox="62 932 559 1071" style="list-style-type: none">• within two hours after a meal• within two hours after a rapid-acting insulin injection or insulin pump bolus, or• during or within two hours after exercise	Fingertip
When you are concerned about the possibility of hypoglycemia	Fingertip

⚠ Do not test on your forearm or palm when:

- ▶ you think your blood glucose is rapidly falling, such as within two hours of exercise or a rapid-acting insulin injection or insulin pump bolus. Testing with a fingertip sample may identify a hypoglycemic (low blood sugar) level sooner than a test with a forearm or palm sample.
- ▶ it has been less than two hours after a meal, a rapid-acting insulin injection or insulin pump bolus, physical exercise, or you think your glucose level is changing rapidly.
- ▶ you are concerned about the possibility of hypoglycemia (insulin reactions), such as when driving a car. This is especially important if you suffer from hypoglycemia unawareness (lack of symptoms to indicate an insulin reaction).

Remember: Consult with your healthcare professional before using your forearm or palm for testing.

Remember, regardless of which sampling site you choose, follow all the instructions in this booklet and your Owner's Booklet.

If you experience symptoms that are not consistent with your test results, and you have followed all instructions for proper testing, call your healthcare professional.

Sampling instructions

These instructions are for getting a sample from your forearm or palm. Your Owner's Booklet contains full instructions for getting a blood sample from your fingertips.

Have these things with you when you test your blood from your forearm or palm:

- ▶ Meter that uses the OneTouch® Ultra® Test Strip
- ▶ OneTouch® Ultra® Test Strips
- ▶ OneTouch® UltraSoft® Blood Sampler
- ▶ Sterile lancet
- ▶ OneTouch® UltraClear® Cap
- ▶ OneTouch® Ultra® Control Solution


The meter, OneTouch® Ultra® Control Solution and OneTouch® Ultra® Test Strips are designed, tested, and proven to work together to produce accurate results. Use no other test strips or control solution with your meter.

Choosing your sample site

Before sampling blood from your forearm or palm, read pages 2 and 3 of this booklet carefully.



Forearm sampling

Choose a fleshy area of the forearm away from bone or visible veins. Select a spot with no hair.

Palm sampling


Choose a fleshy area on the palm below your thumb or pinky. Select a spot with no visible veins and away from deep lines.

Blue cap for fingertip sampling only

OneTouch® UltraClear® Cap is used for forearm and palm sampling only. Replace the blue OneTouch® UltraSoft® Blood Sampler cap with the OneTouch® UltraClear® Cap.

Preparing the OneTouch® UltraSoft® Blood Sampler

⚠ To reduce the chance of infection:

- ▶ Never share a lancet or the OneTouch® UltraSoft® Blood Sampler with anyone.
- ▶ Always use a new, sterile lancet. Lancets are for single use only.
- ▶ Keep the OneTouch® UltraSoft® Sampler clean. Complete cleaning instructions can be found in your Owner's Booklet.

If you do not have a OneTouch® UltraSoft® Blood Sampler, please refer to the instructions that came with your lancing device.

STEP 1

Prepare the sample site

Before you test your blood glucose, wash your hands and the sample site with warm, soapy water. Rinse and dry completely.


Sample in an area with no visible veins, and no hair or lines which may cause your blood sample to smear.

Sometimes there is less blood flow to the forearm than to the fingertips. To help you get a large enough blood droplet, you may gently massage or apply warmth to the site to increase blood flow.

STEP 2

Insert a sterile lancet into the sampler


Remove the sampler cap by twisting it off. Insert the lancet into the holder and push in firmly. Twist the protective disk until it separates from the lancet. Do not twist the lancet.

STEP 3

Install the OneTouch® UltraClear® Cap on the sampler

Twist on the OneTouch® UltraClear® Cap until it is snug.


Replace the regular OneTouch® sampler cap, which is used only for fingertips, with the OneTouch® UltraClear® Cap, which is used for forearm and palm sampling.

STEP 4

Adjust the depth setting


You may have to adjust the sampler to a deeper setting to get a large enough blood droplet from your forearm or palm. Twist the depth adjustment toward the large dots to increase the depth.

STEP 5

Cock the sampler

Slide the ejection/cocking control back until it clicks. If it doesn't click, it may have been cocked when you inserted the lancet.


Getting a blood sample

STEP 1

Lance your forearm or palm

Press and hold the sampler against your forearm or palm for a few seconds, then press the release button.

Keep holding the sampler and cap against your skin until a round droplet of blood forms under the cap. Maintain pressure until the blood droplet is **of at least one microlitre** (● actual size).

STEP 2

Carefully remove the sampler

Gently lift the sampler away from your skin. Do not smear the blood sample.



Important information about getting a blood droplet from your forearm or palm

- ▶ You may need to wait a little longer to get a large enough blood droplet from the forearm or palm. Do not squeeze the site excessively.
- ▶ Do not smear or scrape the droplet with the test strip. If you did not get a big enough droplet, try lancing again in a different area.
- ▶ If the sample droplet runs or spreads due to contact with hair or with a line in your palm, do not try to test that sample. Try lancing again in a smoother area.
- ▶ Remember: You may have to adjust the sampler to a deeper setting to get a large enough blood droplet (● actual size).

If bruising occurs or if you have difficulty getting a sample, consider sampling a fingertip instead. You may want to review the choice of sites with your healthcare professional.

Completing your test

When applying a blood droplet from your forearm or palm to a test strip, keep your palm or forearm steady and bring the meter to the droplet with your other hand.

See the Owner's Booklet for complete instructions for:

- ▶ Applying the blood sample to the test strip.
- ▶ Reading your test result.
- ▶ Error messages your meter may display.
- ▶ Properly disposing of used test strips and lancets.
- ▶ Understanding and managing your results.
- ▶ Caring for your system.

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© 2006 LifeScan, Inc. Milpitas, CA 95035