Bedienungsanleitung/Garantie

Gebruiksaanwijzing/Garantie • Mode d'emploi/Garantie
Instrucciones de servicio/Garantía • Manual de instruções/Garantia
Instruction Manual/Guarantee • Instrukcja obsługi/Gwarancja
Návod k použití/Záruka • A használati utasítás/garancia
Руководство по эксплуатации/Гарантия
Οδηγίες Χειρισμών/Εγγύηση • Kullanma kılavuzu/ Garanti

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Voetmassageapparaat • Masseur pour les pieds
Equipo de masaje para los pies • Aparelho de massajar os pés
Foot massage unit • Urządzenie do masażu stóp
Masážní přístroj pro masáže nohou • Lábmasszázkészülék
Прибор для массажа ног • Συσκευή για το μασάζ των ποδιών
Ayak masaj aleti

FM 2460

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GB

General safety instructions

- **Before commissioning this device please read the instruction manual carefully.**
- Please keep these instructions, the guaranty certificate, the sales receipt and, if possible, the carton with the inner packaging!
- The device is intended exclusively for private and not for commercial use!
- Always remove the plug from the socket whenever the device is not in use, when attaching accessory parts, cleaning the device or whenever a disturbance occurs. Switch off the device beforehand. Pull on the plug, not on the cable.
- In order to protect children from the dangers of electrical appliances, never leave them unsupervised with the device. Consequently, when selecting the location for your device, do so in such a way that children do not have access to the device. Take care to ensure that the cable does not hang down.
- Test the device and the cable regularly for damage. If there is damage of any kind, the device should not be used.
- Do not repair the device yourself, but rather consult an authorised expert.
- For safety reasons a broken or damaged mains lead may only be replaced by an equivalent lead from the manufacturer, our customer service department or a similiary qualified person.
- Keep the device and the cable away from heat, direct sunlight, moisture, sharp edges and suchlike.
- Never use the device unsupervised! Switch off the device whenever you are not using it, even if this is only for a moment.
- Use only original accessories.
- Do not use the device outdoors.
- Under no circumstances must the device be placed in water or other liquid, or come into contact with such. Do not use the device with wet or moist hands.
- Should the device become moist or wet, remove the mains plug from the socket immediately. Do not reach into the water.
- Use the device only for the intended purpose.

Special safety information for this unit

- If the massage feels unpleasant or even begins to hurt, or if you have health problems, please consult your doctor beforehand.
- If you use the unit with liquid, only use water.
- Make sure you are sitting down before you put your feet into the unit.
- Do not use bath oil, gels or other substances.
- Do not use the unit for more than 20 minutes per treatment session.

Operation

GB

1. Put the unit on the floor so that you can sit down comfortably in front of it.
2. Fill up the unit as required with warm or cold water. **Maximum 3.5 litres!**
3. Plug the mains plug into a 230 V, 50 Hz mains socket.
4. Set the operating mode:
 - Whirl pool
 - Whirl pool + heat
 - Whirl pool + heat + massage
 - Intensive massage
 - Off
5. Place your feet in the machine.
6. In order to massage the sole of your foot you can insert three different attachments into the recess in the middle of the machine. These must be inserted crosswise to the angle as seen by the user.
 - The attachment with the large nipples serves for general massage of the feet.
 - The attachment with the small nipples serves as an acupuncture massager.
 - The attachment with the roller is for stimulating and relaxing the sole of the foot.

Switch the unit off after use and pull out the mains plug.

NOTE: The heat produced by the unit is not intended to warm up the water. It is only intended to prevent it from getting cold to quickly.

Cleaning

- Pull the mains plug out of the mains socket each time before you clean the unit.
- Only use a damp cloth to clean the unit.
- Do not use any aggressive or foaming cleaning agents.

This appliance conforms with the CE directives for radio interference suppression and low-voltage safety and has been built to meet current safety requirements. Subject to technical changes without prior notice!

Guarantee

Our appliances are guaranteed for a period of 12 months from the date of purchase (receipt).

Faulty appliances will be repaired or replaced free of charge if these are attributed to faulty material or manufacture.

The guarantee will be invalidated if the appliances are tampered with.

For repairs within the guarantee period, please take the complete appliance with the receipt to the place where it was purchased.

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