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Gebruiksaanwijzing • Mode d'emploi • Instrucciones de servicio
Manual de instruções • Instruction Manual • Instrukcja obsługi
Návod k použití • Руководство по эксплуатации • Kullanma kılavuzu

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Infraroodlamp
Lampe à infrarouges
Lámpara infrarroja
Lâmpada de infravermelhos
Infrared Lamp
Promiennik podczerwieni
Infračervená žárovka
Инфракрасная лампа
Infra kırmızı ışıklı lamba

IRL 2333

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GB

General safety instructions

- **Before commissioning this device please read the instruction manual carefully.**
- Please keep these instructions, the guaranty certificate, the sales receipt and, if possible, the carton with the inner packaging!
- The device is intended exclusively for private and not for commercial use!
- Always remove the plug from the socket whenever the device is not in use, when attaching accessory parts, cleaning the device or whenever a disturbance occurs. Switch off the device beforehand. Pull on the plug, not on the cable.
- In order to protect children from the dangers of electrical appliances, never leave them unsupervised with the device. Consequently, when selecting the location for your device, do so in such a way that children do not have access to the device. Take care to ensure that the cable does not hang down.
- Test the device and the cable regularly for damage. If there is damage of any kind, the device should not be used.
- Do not repair the device yourself, but rather consult an authorised expert.
- Keep the device and the cable away from heat, direct sunlight, moisture, sharp edges and suchlike.
- Never use the device unsupervised! Switch off the device whenever you are not using it, even if this is only for a moment.
- Use only original accessories.
- Do not use the device outdoors.
- Under no circumstances must the device be placed in water or other liquid, or come into contact with such. Do not use the device with wet or moist hands.
- Should the device become moist or wet, remove the mains plug from the socket immediately. Do not reach into the water.
- Use the device only for the intended purpose.

Special Safety Instructions

- **WARNING!** The bulb becomes very hot and requires a long time to cool down. Do not touch. Maintain a distance of 50 cm.
- Only operate the lamp in a standing position on a flat, non-slip surface.
- Please remove all flammable and metal objects from the area to be irradiated.

Using the Lamp

In the case of more serious conditions consult a doctor beforehand.

Operating the Lamp:

1. Place the lamp in position.
2. Insert the plug into a 230 Volt, 50 Hz socket. The bulb will now illuminate.
3. Irradiate the desired area at a distance of 50 cm from the bulb. The heat irradiation will start immediately.
4. The angle of irradiation can be altered by tilting the head of the lamp.

Important Notice

- Irradiation of 15 minutes is sufficient in most cases.
- Irradiation can be repeated as often as desired. However, there should be an interval of several hours between sessions.

Cleaning

- Always remove the **mains plug** before cleaning!
- Clean the **lamp** with a slightly damp cloth.

This appliance conforms with the CE directives for radio interference suppression and low-voltage safety and has been built to meet current safety requirements. Subject to technical changes without prior notice!

Guarantee

Our appliances are guaranteed for a period of 24 months from the date of purchase (receipt).

Faulty appliances will be repaired or replaced free of charge if these are attributed to faulty material or manufacture.

The guarantee will be invalidated if the appliances are tampered with.

For repairs within the guarantee period, please take the complete appliance with the receipt to the place where it was purchased.

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