Terrailon France & Headquarters

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Terrillon®

RELAXING FOOT JACUZZI WITH INFRARED HEAT AND BUBBLES
Instructions and guarantee

EN

BEFORE YOU USE THIS APPLIANCE, IT IS IMPORTANT TO READ THE INSTRUCTIONS AND ADVICE CONTAINED IN THIS BOOKLET CAREFULLY, KEEP THIS BOOKLET IN A SAFE PLACE.

Relaxing foot jacuzzi to soothe aching, tired feet. The product combines the benefits of infrared heat with relaxing massage helping to restore the natural harmony of your entire body. The heating infrared heat penetrates deep into the muscles easing away the aches and stresses of the days. The accessories supplied can be used to target specific reflex zones. The heated base maintains the water at a constant temperature. A relaxing foot massage is a delightful way to end a stressful day.

■ GENERAL RECOMMENDATIONS

The use of any electrical appliance involves respecting certain fundamental rules:

- Unplug the appliance immediately after use, always before cleaning it and before inserting or removing accessories.
- Do not pull the cable or the appliance to remove the plug from the socket.
- If an electrical appliance falls into water, do not try to catch it, but unplug it immediately.
- Do not use the appliance when you are taking a shower or bath.
- Never submerge the appliance in water or any other liquid.
- Do not leave the appliance unattended when it is switched on and do not leave it within reach of children or disabled people.
- This appliance must be used exclusively for the purpose for which it has been designed, and in accordance with the information given in the instruction booklet. Any other use will be considered inappropriate and therefore dangerous.
 The manufacturer cannot be considered liable in the case of damage due to inappropriate or incorrect use.
- Before use, make sure that the appliance is intact and has not been visibly damaged. If you are in any doubt, do not use the appliance and contact the retailer.
- In the case of failure and/or malfunction, switch the appliance off immediately without damaging it. If it needs to be repaired, contact the retailer.
- Keep the appliance away from any source of heat.
- When plugging in or unplugging the appliance and activating the function selection button, make sure that your hands are dry.
- This appliance must be used only with the voltage indicated on the bottom of the product.
- This appliance is not recommended for professional use or treatment.
- After using the appliance for 20 minutes without interruption, it must be left to cool down for at least 15 minutes.

IMPORTANT! This equipment must not be used in a shower room. It must be supplied through a residual differential current device with a breaking current of not more than 30 mA. If the cable is damaged it has to be replaced by the manufacturer, the manufacturer's after-sales service or persons of similar qualifications in order to prevent a hazard.

■ SPECIAL RECOMMENDATIONS

N.B.: The quantity of water must never exceed the maximum level indicated.

- Do not use the appliance in damp places or when using aerosols or administering oxygen.
- Do not use two appliances at the same time under any circumstances (for example: hydromassage foot bath and hair dryer, hydromassage foot bath and electrical razor or electrical toothbrush).
- Do not stand up in the hydromassage foot bath: sit down before you put your feet in the water.
- The appliance cannot be used to heat up cold water. If you wish to have a massage with warm water, fill the appliance with water that is already warm. The temperature will be maintained.
- If you have health problems or have been fitted with a pacemaker, you are recommended to consult your doctor before using the appliance.
- If the massage is unpleasant or if it causes pain both during and after the treatment, it is recommended to consult your doctor. In all cases, you are recommended not to use the appliance if you have any injuries, skin problems, burns or bruises.
- When filling and emptying, the appliance must always be unplugged.

■ INSTRUCTIONS FOR USE

- Make sure that the appliance is unplugged, and then fill the Foot Spa up with warm water.
- Plug in the appliance and select the type of function.

Position

OFF	Off.
Vibration	Gentle vibration massage, Infrared lamp, Maintenance of water temperature.
Vibration & Bubble	Gentle vibration massage, Infrared lamp, Bubbles, Maintenance of water temperature.
Bubble	Bubbles, Maintenance of water temperature

- First, check the temperature of the water, then put your feet in the hydromassage foot bath.
- To empty the appliance, switch it off (position OFF), unplug and empty the water.
- Clean with a damp cloth. Do not use chemical or abrasive agents. Wipe carefully.

■ GUARANTEE

This product is guaranteed for 2 years against any defects in materials or manufacturing. During that period, any such defects will be repaired free of charge (proof of purchase must be shown in the case of a claim under warranty).

The warranty does not cover damage resulting from accidents, improper use or negligence. In the case of a claim, first contact the shop where you purchased this product.

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