Français

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GB English Instructions and warranty

DIGITAL THERMOMETER

Before use, read the instructions and warnings given in this handbook carefully and keep the handbook in a safe place. Digital thermometer with a large LCD display for easy reading.

Fast and accurate reading with a sound signal when the temperature reading has been taken. Stores the last read temperature, to evaluate any variations. Extremely easy to use, it is kept in a practical case to avoid risk of breaking. The probes can be washed, and it turns off automatically.

■ GENERAL WARRINGTON

Self-measuring means control, not diagnosis and treatment, any unusual readings must be discussed with your doctor. Under no circumstances must you change the medicine doses that have been prescribed by your doctor.

2 Do not drop, bend, drop or break the thermometer.

3 Keep the thermometer probe clean before and after each use.

4 Do not expose to sprays or direct sources of heat. Keep in a dry clean place.

5 Do not immerse the thermometer in water or other liquids, only the probe is waterproof.

6 The thermometer reads only the rectal, oral or oral temperature. Do not use it for reading armpit temperatures.

■ INSTRUCTIONS FOR USE

1 Press the ON/OFF key. For a few seconds all the function symbols light up on the display. After the thermometer displays 36,5°C, when the °C symbol starts flashing, the thermometer is ready for use.

2 Place the probe as follows: IN THE MOUTH: place the probe under the tongue and keep the mouth closed for about 2-3 minutes, until the temperature stabilises. The standard temperature for this type of reading is 36,7°C. UNDER THE ARM: place the probe under the armpit and keep the arm well pressed against the chest for about 2-3 minutes, until the temperature stabilises. The standard temperature for this type of reading is 36,5°C. IN THE RECTUM: lubricate the probe with a water-based lubricant. Insert the probe into the rectum for about 1 min., until the temperature stabilises (about 2-3 minutes). The standard temperature for this type of reading is 37,0°C.

3 When the thermometer gives off the sound signal, remove it. The reading is terminated and shown on the display.

4 The reading is terminated automatically after 10 minutes without use (in the meantime the reading is still displayed) otherwise you can turn it off by pressing the ON/OFF button.

If the temperature is above 44°C, the display shows HC; if it is below 32°C, the temperature shows LC.

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CE 0434 (Dir. 93/42/EEC)

■ MEMORY FUNCTION

The thermometer stores the last reading in its memory, to evaluate any variations.

To recall the last reading, press the ON/OFF button for a few seconds until the reading appears on the display, in the bottom right corner of the display, "M" appears.

Release the ON/OFF button to continue with measuring the actual temperature.

■ MAINTENANCE

Clean the thermometer probe with a cloth dampened with a 0.02% hypochlorite sodium solution or with alcohol or with any type of topical disinfectant before and/or after each use. The rest of the thermometer must be cleaned with a dry cloth. Keep the thermometer safe in its special case.

■ TECHNICAL DATA

Measurement range: 32°-44°C
Tolerance: +/- 0,1°C 35-39°C, otherwise +/- 0,2°C
Accuracy: +/- 0,1°C

Working conditions: 0°C to 50°C, Humidity: 10-83%

Conservation conditions: -20°C to 60°C, Humidity: 10-83%

Battery autonomy: approx. 100 hours continuous use

Dimensions: 126,5 mm x 18 mm x 10 mm

Weight: 7 g (with battery)

■ WARRANTY

This appliance is guaranteed for a period of 2 years from the date of purchase as shown by the receipt issued by the supplier and by the manufacturer or by the distributor for whom which must be kept attached to this warranty. The warranty only covers damage in materials and workmanship and is not applicable if damage results from an accident, incorrect use, negligence or improper use of the product. Only those accessories supplied with any other accessories may cause the warranty to be void. If there is a fault, contact the dealer: do NOT send equipment directly to the manufacturer. **No repairs or returns will be accepted.** The company of the manufacturer does not accept responsibility for any damage to persons or pets or damage to property, whether direct or indirect, resulting from non-compliance with all provisions indicated in the relevant instruction book, especially concerning warnings regarding installation, use and maintenance of the appliance.

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