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ELECTRONIC BATHROOM SCALE - BODY FAT ANALYSER

Dear Customer,

Thank you for purchasing this product. We hope you will find it satisfactory. Please read these instructions carefully to take full advantage of your product.

READ BEFORE USE

- Do not use the scale on subjects who have a pacemaker. It may cause the device to malfunction. When in doubt, please consult your physician.
- For pregnant women and people who have body implants such as artificial limbs, contraceptive devices, metal plates or screws, the results given by the Body Fat Analyser may be inaccurate. However there is no health risk associated with using the Body Fat Analyser.
- Do not disassemble the scale as incorrect handling may cause injury.
- The product is not intended for professional usage in hospitals or medical facilities; it is intended for home/consumer use only.
- Clean after use with a lightly dampened cloth without allowing any water to get into the device. Do not use solvents or immerse the product in water.
- Avoid excessive impact or vibration to the unit, such as dropping the product onto the floor.
- The scale must be reset each time it has been moved.
- Pause after each weighing.

PRODUCT SPECIFICATION

function keys

3-BUTTON MODEL		4-BUTTON MODEL	
	<ul style="list-style-type: none"> • Turns the unit on and off • Confirmation button for Body Fat function, toggles from Height, Age, Gender and Athlete/Normal for Body Fat analysis • Recalls the 4 parameters (Height, Age, Gender & Athlete/Normal) from 4 memories 	ON	<ul style="list-style-type: none"> • Turns the unit on and off • Recalls the 4 parameters (Height, Age, Gender & Athlete/Normal) from memories
			<ul style="list-style-type: none"> • Confirmation button for Body Fat function, toggles from Height, Age, Gender and Athlete/Normal for Body Fat analysis
	<ul style="list-style-type: none"> • Changes value of height and age, toggles between male, female and between Athlete/Normal for the Body Fat function • Selects memory in Memory mode & Recall mode • Clears contents of the user memory in Recall mode 		<ul style="list-style-type: none"> • Changes value of height and age, toggles between male, female and between Athlete/Normal for the Body Fat function • Selects memory in Memory mode & Recall mode • Clears contents of the user memory in Recall mode

Screen Display

	Height		Body Water Percentage
	Weight		Bone Mass
	Body Fat Percentage		Muscle Mass
N/ 	Normal mode	OC	Corresponds to the daily calories intake you need to ingest to achieve your ideal weight which depends on your age, size, gender and level of physical activity
A/ 	Athlete Mode. An athlete is defined as a person who consistently trains a minimum of three times per week for two hours each time		
PC	Corresponds to the daily calories intake you need to ingest to maintain your actual weight		Corresponds to the activity level (1 to 5, 5 being the highest activity level)

ELECTRONIC BATHROOM SCALE - BODY FAT ANALYSER

PRODUCT USAGE

Initialisation

- Insert the battery(s) in the battery compartment, making sure that the polarities are correct, and/or pull the plastic tape to initialise the batteries for the first time. Remove the battery(s) if the product is not in use for an extended period. The battery type is specified inside the battery compartment.
- The scale needs to be initialised at the first time of use or after battery replacement. Place the scale on a hard, flat surface. Press briefly and sharply with your foot in the middle of the platform. The display reads OOOO or - - - -, then the scale automatically turns off. Your scale is now ready for use.
- If your scale has a conversion switch, you can change the measuring units from metric to imperial or vice versa (i.e. KG/CM - ST/IN - LB/IN) by simply adjusting the conversion switch located on the under-side of the scale.

Note: it is recommended to take weight and body fat measurements at a constant time period daily for example: from 6 am to 7 am daily.

Operation

Please follow the illustrated directions at the beginning of the instruction manual.

A : Weight measurement only

B : Weight measurement with BFA function: measure body fat, body water percentage, muscle and bone mass

C : Memory : save (C.1), quick recall (C.2)

Warning indicators

Err0	Incorrect initialisation – The scale must be reset	Err4	Percentage of body fat is outside the 5% to 50% measurement range: begin again to check
Err1	Unstable weight: begin again	Err5	Percentage of body water is outside the measurement range: clean the bottom of feet with damp cloth and step on the scale again to check
Err2	Overload: remove weight immediately	Lo	Battery failure: replace it
Err3	Body impedance exceeds the limit (100 ohm - 1kohm): begin again. For certain skin types an Err3 reading may occasionally appear. In this case you may need to dampen your feet to improve the contact with the electrode		Full memory. Press (▼/▲) to clear

PROTECT THE ENVIRONMENT

Dispose of used batteries in a designed container so that they can be collected and recycled. Do not mix different types of batteries. Do not mix new and used batteries.

After life, dispose of this product at a designated waste recycling point.

GUARANTEE

This product is guaranteed against defects in materials or manufacturing. During the guarantee period, any such defects will be repaired free of charge (proof of purchase must be shown in the case of a complaint under guarantee). The guarantee does not cover damage resulting from improper use or negligence. In the event of a claim, first contact the shop where you purchased the product.

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