Haute technologie
High technology

Nouveaux courants
New currents

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FACE
4 électrodes
rechargeables

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- Design patent pending
- US Patent n° 6,122,554

FEATURES COMMON TO

SPORT-ELEC Evolution Body and SPORT ELEC Evolution Body&Face

General

SPORT-ELEC « Evolution » is a passive gymnastic apparatus using electrical stimulation through the "contraction-release" effect on muscles that you want to work. Our bodies contain more than 600 muscles. The shape and harmony of our bodies depends on these muscles. **SPORT-ELEC « Evolution »** allows you to activate a set of muscles, without any effort from you. Its daily use will help you to not only **refine, remodel your silhouette, maintain and tone your musculature, but also relax and enjoy the beneficial effects of a massage.**

SPORT-ELEC « Evolution » diffuses a very high quality current - the two-phase current. This is a current tested and recognised to guarantee optimal comfort and safety in use.

SPORT-ELEC « Evolution » contains 5 programmes for body and face resulting from **SPORT-ELEC** research. The frequencies and pulse widths of each programme have been defined to correspond perfectly with each person's needs.

New design: New harmonious shape, compactness, handiness, practicality - all combine to form the top advantages of **SPORT-ELEC « Evolution »**.

Leading edge technology: Extreme care is given to the current's quality and the definition of programmes.

Easy to use: Fitted with a micro-controller, a timer (automatic stop at the end of each programme) and an automatic programme selector: **SPORT-ELEC « Evolution »** is truly accessible to all.

Choose your programme. **SPORT-ELEC « Evolution »** looks after you.

NEW !

Thanks to a completely new system - the performance index - **SPORT-ELEC « Evolution »** allows you to assess your performances and the effectiveness of your session. (see "performance control" function)
By choosing **SPORT ELEC**, you are trusting yourself to a major brand which is recognised in the electro-beauty field and complies with all safety and quality standards required in Europe.

SPORT ELEC EVOLUTION Body Presentation

See photo page 2

- 1- LCD screen: programme identification, performance indicator
- 2- "Performance control" Index (performance assessment figure)
- 3- "programme selection" button
- 4- Power-on button for Dark Blue channel (channel A)
- 5- Power-on button for Light Blue channel (channel B)
- 6- On-off switch
- 7- 80 cm long flexible wires
- 8- 50 mm diameter body electrodes
- 13- Wire guide
- 14- Plug

Presentation of SPORT-ELEC EVOLUTION Body and Face

See photo page 2

- 1- LCD screen: programme identification, performance indicator
- 2- "Performance Control" index (performance assessment figure)
- 3- "Programme selection" button
- 4- Power-on button for Dark Blue channel (channel A)
- 5- Power-on button for Light Blue channel (channel B)
- 6- On-off switch
- 7- 80cm long flexible wires
- 8- 50mm diameter body electrodes
- 9- Built-in face connection (also available on **EVOLUTION Body**)
- 10- 4 60 x 20 mm rechargeable face electrodes with cosmetic (**EVOLUTION serum**) Patented system
- 11- 1 30 ml/1 Floz tube **EVOLUTION Serum**
- 12- 1 30 ml/1 Floz tube **EVOLUTION Gel**
- 13- Wire guide
- 14- Plug
- 15- Hook for holding face electrode wires.

SPORT-ELEC « Evolution » is a registered design. It is presented in the form of a case in shiny ABS.
Two groups of 3 wires (2 dark blue + 1 grey, 2 light blue + 1 grey) 80 cm in length connect the **SPORT-ELEC « Evolution »** appliance to the body electrodes (50mm in diameter).

CONTRA-INDICATIONS :

- Do not use **EVOLUTION :**
- On the cardiogenic area
- If you have been fitted with a pacemaker
- If you are pregnant

SAFETY PRECAUTIONS

- * Do not use body electrodes on the face
- * Do not leave within the reach of children
- * Ask for advice from your doctor :
 - In case of electronic or metal implants in general (pins, etc..)
 - In case of dermatitis (wounds, etc..)
- * The appliance is not recommended for therapeutic use (seek advice from your doctor).

WARNING

Do not use your **SPORT-ELEC « Evolution »** in a damp atmosphere (bathroom, sauna) or close to an electrical appliance connected to the mains and the ground, or even piping connected to the ground. A patient's simultaneous connection to high-frequency surgical equipment may cause burns at the contact points of the stimulator's electrodes and the stimulator may possibly be damaged. Operation in the immediate vicinity of a short-wave therapy apparatus may cause instabilities in the power output from the stimulator.

IMPORTANT : IDEAL TO USE A MONTH AND A HALF AFTER GIVING BIRTH to help you improve muscle tone.

USING THE SPORT-ELEC EVOLUTION for the body

- 1° Check that the switch for the **SPORT-ELEC « Evolution »** is on "0" and that the power-on buttons are at minimum (pushed to the left). Open the "battery" compartment and insert 2 9 V batteries into the compartment and then close the case. Comply with + and - indications.

(See diagram A Page 4)

1. Gripper clip
2. Battery compartment
3. 9 volt batteries (not provided)
4. Detachable cover
Battery life: 1 month if used for 20 mn/day. If not used for long periods, remove batteries from case.

- 2° Take the 6 electrodes and carefully remove the holder which protects each of them, moisten them one by one, spreading a few drops of water on the adhesive gel with the tips of your fingers and wait for 2 minutes.

(see diagram B page 4)

- 3° The electrodes must be placed symmetrically, for example: on the stomach - to the left, 2 dark blue and 1 grey connections, to the right, 2 light blue and 1 grey connections. Place the electrodes on previously washed and dried skin, always respecting a "triangular" arrangement.

(see diagram C page 4):

G = grey, A = dark blue, B = light blue

Position the electrodes on your body as shown in diagrams D-E-F-G-H-I-J-K page 4:

G = grey, A = dark blue, B = light blue

D - Pectoral muscles	10 - Middle trapezius
1 - Major pectoral	11 - Lower trapezius
α - Do not place electrodes in the heart region	I - Back and buttocks
E - Shoulders	12 - Main back muscles
2 - Middle deltoid	13 - Main gluteal muscles
3 - Front	J - Legs
4 - Back	14 - Crural biceps
F - Abdominal muscles	15 - Half tendinous
5 - Major oblique	16 - Outer gemellus
6 - Major longitudinal	17 - Inner gemellus
G - Biceps	K - Thighs
7 - Biceps - longitudinal	18 - Major adductors
8 - Biceps - transversal	19 - Vastus internus
H - Trapezii	20 - Rear longitudinal
9 - Upper trapezius	21 - Vastus externus

- 4° Attach the **SPORT-ELEC « Evolution »** to the edge of an item of clothing or underwear with its clip.

- 5° Take the wire guide and bring it up to the plugs so as to separate the dark blue + grey group from the light blue + grey group. Connect the electrodes to the appliance as shown in the following pages.

Specific features of the SPORT-ELEC EVOLUTION Body&Face

Thanks to its face electrodes + serum and gel, **SPORT ELEC Evolution Body & Face** makes it possible to use the "Face Toning" programme. The aim is to maintain the firmness of the oval of your face, attenuate the effects of tiredness, by working on the areas around your eyes and forehead, to tauten the lower face by preventing the "cheek-pouch" effect. This programme offers a full face treatment which is hydrating, revitalising and tonic.

SPORT ELEC EVOLUTION Body&Face is the only apparatus to offer the combination of electrical + cosmetic stimulation of muscles thanks to its patented system of electrodes.

USING THE FACE PROGRAMME

Impregnate the face electrodes with a dozen drops of **SERUM** and then remove the electrodes from their plastic holder. (see diagrams p 7-a.)

Slide the wires of the face electrodes into the hook and fix it to the edge of a garment.

Place the electrodes (see photos M.N p7) on your face on previously washed and dried skin, plug your face electrodes into the connector provided for that purpose (see photos c p7)

Select the **FACE TONING** programme, gently increase the power of each channel (repeat the same procedure for use detailed on page 16), leave it to act. You can work both sides of your face at the same time (2 electrodes on one side, 2 electrodes on the other).

USER TIPS : «FACE TONING»

Gently turn the dark blue button for channel A from 1 to 8, you will begin to feel contractions, increase power as necessary. The contractions must be firm but never painful. Proceed in the same way for the light blue power-on button for channel B.

IMPORTANT : DO NOT USE ANY OTHER COSMETIC PRODUCT THAN THAT SPECIFIED FOR FACE ELECTRODES.

COMPOSITION OF ACTIVE INGREDIENTS IN SERUM AND GEL (see ingredient list p.7)

SPORT-ELEC EVOLUTION SERUM :

ACTIVE INGREDIENTS:

- **Tensor action:** wheat proteins. Smooths the skin's micro-relief, helps to procure the effect of "healthy glow".
- **Moisturising action:** Sodium PCA. Helps to maintain the hydration of the outer layers of the epidermis.
- **Revitalising action:** Seaweed extracts. Rich in mineral salts and trace elements.
- **Tonic action:** yeast extracts.

SPORT-ELEC EVOLUTION GEL:

applied once the session has finished, has an immediate tensor action, and vitamin-producing action, gives you back a truly radiant skin.

ACTIVE INGREDIENTS:

- **Immediate tensor action:** Wheat proteins. Smooths the facial features and immediately restores radiance to your face.
- **Vitamin-producing action:** Fruity lemon-orange complex. Rich in vitamin C.

The face electrodes must be plugged into the face connector provided for that purpose.

For best results, we recommend one 40 minute session every day for 4 weeks. Enough of each of the products is supplied for 30 treatment sessions. Of course, the machine can be used for more than one 30-day course with the Evolution Recharge Kit which contains 2 electrodes, 1 bottle of Serum and 1 bottle of Tensor Gel. For further information, contact your local SPORT-ELEC supplier or phone to : +33-(0)2-32-96-50-50.

 **NB: The RELIEVE RELAXE MUSCULAR FATIGUE programme and the FACE TONING programme are represented by the same pictogram.**

N.B. DO NOT USE THE BODY ELECTRODES AT THE SAME TIME AS THE FACE ELECTRODES.

"PERFORMANCE CONTROL" FUNCTION EVOLUTION has a performance index which allows you to assess your results.

Taking into account the real duration and the level of power selected, **SPORT-ELEC « Evolution »** displays a figure which will calculate the level of effectiveness for your session. This figure appears on the bar-graph shown at the bottom of the screen, at the end of the session. Example: The bar-graph displays 6. The following table helps you to assess the level of effectiveness and performance corresponding to this figure.

- 1: Cool
- 2: Calm
- 3: Amateur
- 4: Sport
- 5: Super Sport
- 6: Pro
- 7: Super Pro
- 8 : Champion
- 9 : Super Champion

7° Once the programme has ended, return the power-on buttons to minimum and move the switch back to the "0" position. Gently disconnect and remove the electrodes, moisten them slightly with a few drops of water and put them carefully back on their holder. Electrodes are stored in their original packaging in a cool place (vegetable compartment of a fridge).

8° USER TIPS

- It is possible to restart the cycle of one of the programmes on another muscular zone without however exceeding the 40 min of use on the same part of the body.
- If you wish to stop a programme in progress, all you have to do is put the 2 buttons back to "0" and repeat the procedure for starting up the appliance. Then select a new programme and repeat the procedure for starting up the appliance.
- For optimal use, take inspiration from the diagrams reproduced in this brochure.
- The case and wires of **SPORT ELEC** may be cleaned with a blotter moistened with a disinfectant solution (such as Mercryl) diluted x 10 with water)
- Parts in contact, such as electrodes, may be rehydrated before and after use with a disinfecting agent (e.g. Hèxomédine)
- The electrodes must be repositioned on their holder and stored in a cool place between +5°C and +10°C.
- If not using **SPORT ELEC** for any length of time, please remove batteries from case.

The 6 electrodes have a life of approximately 40 applications, don't hesitate to order them from us directly at the company's head office.

SPORT-ELEC S.A.
B.P. 35 - 27520 BOURGTHEROULDE
TEL : + 33 (0)2.32.96.50.50
 or contact your local supplier

6° START OF SESSION/PROGRAMME SELECTION

a / Make sure that the two buttons for the Dark Blue (channel A) and Light Blue (channel B) channels are at minimum (safe position for starting up the appliance).

Move switch to "I".

b / CHOOSE YOUR PROGRAMME:

All you have to do is press on the "programme selection" button placed under the screen, and the programmes will scroll:

FEEL-GOOD / WARM-UP MASSAGE

20 Min for warming up your muscles before sport, or relaxing after an effort.

FIRM-UP STIMULATION

This 40 min programme offers the possibility of effectively firming up your body after a period of inactivity, and maintaining your body's tonus.

TARGETED BODY ACTION

Harmonises your silhouette. This 20 min programme is ideal for specifically working on an area of your body that you particularly want to improve.

MUSCULAR FATIGUE: relieves, relaxes

FACE TONING: (see "SPORT-ELEC Evolution Body & Face")

c / ADJUSTING INTENSITY:

Gently turn the dark blue button for channel A from left to right, you will begin to feel contractions, increase power as necessary. The contractions must be firm but never painful.

Proceed in the same way for the light blue power-on button for channel B.

If you feel an unpleasant sensation around the electrodes: bring the power-on buttons back to minimum, move the electrode(s) and gradually increase the power. Once the adjustments have been optimised in this way, leave your appliance to work for the full duration of the programme.

With the "tri-electrode" system, it is possible that the sensation obtained will not be the same across the 3 electrodes. You may have a sensation from only one electrode. In this case, move the electrode attached to the grey connector so as to find the best possible contraction with the minimum sensation. The quality of the contraction remains the same.

Recommendations :

- The electrodes and related equipment are environment-friendly and can be placed in household waste bins for normal disposal.

- The used batteries should be disposed of in existing collecting systems or special containers.

DESCRIPTION OF PROGRAMMES

FEEL GOOD / WARM-UP MASSAGE **TOTAL TIME : 22'**
REPETITION : 20x

Phase	Frequency	Pulse Width	Duration
Phase1	1Hz	2.0µs	4"
Phase2	2Hz	2.0µs	4"
Phase3	3Hz	2.0µs	4"
Phase4	4Hz	2.0µs	4"
Phase5	5Hz	2.0µs	4"
Phase6	6Hz	2.0µs	4"
Phase7	7Hz	2.0µs	4"
Phase8	8Hz	2.0µs	4"
Phase9	9Hz	2.0µs	4"
Phase10	10Hz	2.0µs	4"
Phase11	11Hz	2.0µs	4"
Phase12	12Hz	2.0µs	4"
Phase13	13Hz	2.0µs	4"
Phase14	14Hz	2.0µs	4"
Phase15	15Hz	2.0µs	4"
Phase16	16Hz	2.0µs	4"

FIRM-UP STIMULATION **TOTAL TIME : 30'**
REPETITION : 26x

Phase	Frequency	Pulse Width	Duration
Phase1	35Hz	2.2µs	10"
Phase2	8Hz	2.2µs	6"
Phase3	45Hz	2.8µs	10"
Phase4	8Hz	2.2µs	6"
Phase5	55Hz	2.8µs	10"
Phase6	8Hz	2.2µs	6"
Phase7	65Hz	2.8µs	10"
Phase8	8Hz	2.2µs	6"

R = Recovery period

TARGETED BODY ACTION **TOTAL TIME : 33'**
REPETITION : 30x

Phase	Frequency	Pulse Width	Duration
Phase1	75Hz	2.8µs	4"
Phase2	3Hz	2.8µs	12"
Phase3	85Hz	2.8µs	4"
Phase4	3Hz	2.8µs	12"
Phase5	95Hz	2.8µs	4"
Phase6	3Hz	2.8µs	12"
Phase7	105Hz	2.8µs	4"
Phase8	3Hz	2.8µs	12"

R = Recovery period

MUSCULAR FATIGUE : relieves, relaxes **TOTAL TIME : 20'**
REPETITION : 16x

Phase	Frequency	Pulse Width	Duration
Phase1	40Hz	1.2µs	1.2"
Phase2	41Hz	1.2µs	1.2"
Phase3	42Hz	1.2µs	1.2"
Phase4	43Hz	1.2µs	1.2"
Phase5	44Hz	1.2µs	1.2"
Phase6	96Hz	1.2µs	1.2"
Phase7	97Hz	1.2µs	1.2"
Phase8	98Hz	1.2µs	1.2"
Phase9	99Hz	1.2µs	1.2"
Phase10	100Hz	1.2µs	1.2"

FACE TONING (SPORT-ELEC Evolution Body & Face) **TOTAL TIME : 20'**
REPETITION : 16x

Phase	Frequency	Pulse Width	Duration
Phase1	40Hz	1.2µs	1.2"
Phase2	41Hz	1.2µs	1.2"
Phase3	42Hz	1.2µs	1.2"
Phase4	43Hz	1.2µs	1.2"
Phase5	44Hz	1.2µs	1.2"
Phase6	96Hz	1.2µs	1.2"
Phase7	97Hz	1.2µs	1.2"
Phase8	98Hz	1.2µs	1.2"
Phase9	99Hz	1.2µs	1.2"
Phase10	100Hz	1.2µs	1.2"

PHASE : Time of diffusion of an electrical current.

REPETITION : Number of phases which are repeated inside each programme. The figures indicating the number of repetitions are rounded up to the

TOTAL TIME : Length of time for each programme.

CLINICAL TESTS

Clinical tests were conducted in France by an accredited institution.

Subjects : 30 adult, female volunteers of between 26 and 64 years of age for evaluation of Sport-Elec's moisturising action following a single session. 19 adult, female subjects of between 35 and 64 years of age for evaluation of Sport-Elec's firming, tensor, toning and moisturising actions after daily use of the apparatus for a period of 3 weeks. 10 adult, female volunteers of between 37 and 64 years of age for evaluation of Sport-Elec's action against wrinkles (folded skin). Results recorded after 30 minutes daily use of Sport-Elec for a period of 3 weeks (under medical supervision):

- A clear effect on the moisture content of the skin (the outermost layers of the epidermis)
- A reduction of the order of 23 % in the number of deep wrinkles*.
- A reduction of the order of 30 % in the actual depth of the wrinkles*.

* The figures are representative but the product's efficacy varied between different subjects (depending on skin-type, age, etc.). Evaluation of the product's acceptability by volunteers who had used it for 4 weeks: Acceptability: Good to excellent 100%.

9° SAFETY CHECK

Every 24 months, the following points must be checked:

- 1- Check whether there has been superficial deterioration of the material
- 2- Check that the appliance is functioning correctly:
 - * legibility of writing on the case
 - * switch, contactor
- 3- Check the reliability of accessories:
 - * wire from electrodes in good condition
 - * electrodes in good condition (no cracks, no defects in material)
- 4- Measure electrical values:
 - * for a charge of 500 ohms, the peak voltage at outlet must not exceed 100 V.

10° ENVIRONMENTAL CONDITIONS FOR TRANSPORT AND STORAGE:

- Ambient temperature between +10°C and + 40°C
- Relative humidity between 30% and 75%
- Pressure 700 hPa

Technical characteristics of the SPORT ELEC Evolution

- Two-phase current
- Battery power supply : 2 x 9 v (not provided)
- Battery type : 6LR61 – 6F22
- Pulse frequency : 1 hz to 100 Hz
- Pulse width : from 120 µs to 280 µs
- Adjustable intensity, 2 independent generators
- DIMENSIONS : 75x145x40 mm
- WEIGHT : 190 g

80 Vcc/500 Ω

Standard AF classification

N.B. read operating instructions carefully

European directive (89/336/CEE) CEM (76/768/CEE) Cosmetic

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