
ES - MANUAL DE INSTRUCCIONES
PT - MANUAL DE INSTRUÇÕES
EN - INSTRUCTIONS FOR USE
FR - MANUEL D'UTILISATION
EL - ΕΝΤΥΠΟ ΟΔΗΓΙΩΝ

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1. INTRODUCTION

This appliance works on the principal of measuring the body's electronic impedance, which is referred to using the term bioelectric, i.e. "bioelectric impedance analysis" (BIA), a factor relating to the percentage of body fat and body water, while taking other physical factors into consideration (age, sex and weight). This weighing scale can also measure bone mass, muscle mass and the number of kilocalories required during the day.

Before using the scales for the first time, read this manual carefully to guarantee its correct functioning and extend the life of the appliance. Keep the instructions for later reference.

2. IMPORTANT INFORMATION

- The readings for your body fat and water percentage are guidance measurements only. This appliance is not recommended for therapeutic use. Please consult your doctor for advice on any medication or diet you should follow in order to reach a normal weight.
- The use of these scales is not recommended for persons using pacemakers or other therapeutic implants with electronic components. It could be dangerous and could reduce the performance of the implant.
- To avoid any interference, do not use this appliance near devices with a strong electromagnetic field. Always keep well away from any such devices when you are using the scales, or simply switch off any devices emitting strong electromagnetic fields before using it.

3. APPLIANCE DESCRIPTION (FIG. 1)

- Display.
- \blacktriangle and \blacktriangledown : Data selection buttons.
- Set button.
- Weighing surface.
- Battery compartment.

4. CHARACTERISTICS / SPECIFICATIONS

- The scales are supplied programmed with the following data:

Parameter	Default programming	Range
Memory	1	1-12
Sex	Male	Male/Female
Height	170 cm	80 cm - 220 cm
Age	30	6 -100

- Maximum weighing capacity: 150 kg / 330 lb.
- Minimum: 3 kg / 6.6 lb.
- Graduation: 0.1 kg / 0.2 lb.
- Weight tolerance: $\pm 1\%$
- Start-up and automatic switch-off.
- Low battery indicator ("LO")
- Overload indicator (O-Ld)
- Power source: Two 3V lithium batteries (CR 2430)

5. SCREEN DISPLAY (FIG. 2)

- Male
- Female
- Memory
- Height
- Age
- Weight unit
- Body water percentage
- Status line
- Body fat percentage
- Bone mass
- Muscle mass
- Kilocalories

6. ADJUSTMENT

1. Insert the two batteries in the battery compartment (5) on the underside of the scales.
2. If you only wish to know your weight, press one of the upper corners of the scale and select the weight unit (kg or lb) using the Set button (3).
3. If you wish to use the scales in calculation and analysis mode, follow the instructions in point 7.
4. Place the scales on a hard, flat surface (do not place them on a carpet). To guarantee maximum accuracy and repeatability of the readings, always weigh yourself on the scales in the same place.

kg

- III. The next parameter to be entered is **sex** (a-b). When the symbol starts to flash, use the **▲** and **▼** buttons (3) to select male or female. Then press "Set" (4) to confirm the data.


- IV. Another number will then begin to flash, indicating the person's **height** (d). Use the **▲** and **▼** (2) buttons to select the height. Then press "Set" (3) to confirm the data.

176 cm

- V. The last data item is **age** (e). Enter your age using the **▲** and **▼** buttons (2). Then press "Set" (3) to confirm the data.

AGE 25


To make sure all the parameters have been correctly programmed, press "Set" (3) consecutively when you finish programming.

- I. To switch on the scales, press "Set" (3). When the number beside the icon begins to flash, you can select the desired memory setting (c) by moving up and down with the **▲** and **▼** buttons (2). Then press "Set" (3) to confirm the data.
- II. Then select the **weight unit** (f). To change from one weight unit to another, use the **▲** and **▼** buttons (2), and then confirm with the Set button (3).

When you wish to weigh yourself:

1. If the scales are switched off, press "Set" (3). The following will appear:

You do not need to enter your data each time you switch on the scales in calculation and analysis mode. Select the memory position in which you have entered your personal data using ▲ and ▼ (2), and wait for the number to stop flashing.

The scales are ready when the following appears on the display:

2. Stand on the scales carefully, barefoot and with dry feet. Dry the weighing surface for better conduction. Stand on the scales with your feet apart and in line.
3. As soon as you stand on the scales the digits on the display will increase. When the digits stop increasing and begin to flash, the scales are showing your real weight.
4. Stand still on the scales while they calculate your weight values. After a few seconds, the scales will show the number of kilocalories or amount of energy you need per day (l), your bone mass (j), muscle mass (k), body water percentage (g) and body fat percentage (i) and, lastly, your body status (h) (1 mark → underweight, 2 marks → healthy, 3 marks → overweight, or 4 marks → very overweight).


5. The scales will switch off automatically after 20 seconds.

The status line at the bottom of the display, with segments shown in vertical bars (h). This line gives an indication of body fat level, showing how your level compares to the statistical standard, and can be taken as a weight status reference.

NOTE: Whenever you press "Set" (3) and switch the scales on again, you can follow steps I – V again to check or programme all your personal and previous measurement reports before taking new measurements.

WEIGHT MODE PROGRAMMING

Step carefully onto the scales (do not jump onto them).

As soon as you stand on the scales the digits on the display will count up and stabilise at your weight. After a few seconds the display will automatically switch off.

NOTE: For correct weighing, weigh yourself without clothes and without shoes. Always weigh yourself as far as possible in the same circumstances (in the morning before breakfast for example) and at the same time of day.

9. BATTERY LOW INDICATOR

When the battery needs replacing, the following message will appear on the display:

LO

When this happens, follow the adjustment instructions and replace the battery.

Remember that batteries are not household waste. Spent batteries must be disposed of in the special containers provided.

10. REFERENCE INFORMATION

Body fat evaluation is based on the standard method of measurement. It is calculated using a correlated formula developed on the basis of extensive scientific studies. The water content in the muscle without fat affects the result of the measurement. Your particular habits and lifestyle affect the body water level balance, as do other factors such as the following:

- The reading may be lower if you have just taken a bath.
- The readings may be higher after eating and drinking.
- For women, the reading may change during your monthly period.
- Body water may be lost due to illness or physical activity causing effusive sweating and dehydration.
- The following table shows physical status in relation to sex and age.

FEMALE

Age	Body fat				Body water percentage	Muscle mass	Calories per day
	Underweight	Healthy	Overweight	Obese			
6-12	<7	7-19	19-25	>25	>58	>36	1050-2300
13-19	<20	20-31	31-37	>37	>56	>34	
20-39	<21	21-32	32-38	>38	>54	>33	
40-59	<23	23-34	34-40	>40	>53	>32	1125-2700
60-79	<25	25-36	36-42	>42	>52	>30	

MALE

Age	Body fat				Body water percentage	Muscle mass	Calories per day
	Underweight	Healthy	Overweight	Obese			
6-12	<7	7-18	18-25	>25	>60	>40	1175-2800
13-19	<8	8-19	19-25	>25	>58	>38	
20-39	<9	9-20	20-26	>26	>56	>37	
40-59	<11	11-22	22-28	>28	>55	>36	1450-3300
60-79	<13	13-24	24-30	>30	>54	>34	

11. ERROR AND OTHER MESSAGES ON THE DISPLAY

1. The display may occasionally show the error message "E". If this happens, either switch off the scales for a moment or remove the battery for complete reprogramming. Then switch on the scales again for them to resume functioning.
2. "LO": Low battery indicator.
3. "O-Ld": Overload indicator.

12. MAINTENANCE, CLEANING AND CORRECT USE

- The scales may be slippery when wet. Make sure they are dry before using them.

- Clean the scales with a dry cloth and make sure no water enters them. Do not use strong cleaning chemicals such as alcohol or benzene.
- Keep the scales in a cool dry place, rather than in a damp or hot place.
- Keep the scales in horizontal position rather than vertical, as the batteries may run down.
- Do not attempt to open the scales to repair them, as this will nullify the guarantee.
- Do not drop the scales on the floor or jump onto them when you use them. Avoid any knocks and do not shake the scales, as they are a precision appliance. Their structure is delicate and sensitive to any external impact.
- If the scales do not work when you switch them on, make sure a battery has been inserted and that it has not run down. Insert or replace the battery and consult your supplier if necessary.

13. READ THE FOLLOWING NOTES CAREFULLY FOR A CORRECT UNDERSTANDING OF THE APPLICATION AND THEORY OF THE APPLIANCE

1. The body fat and body water content measured using the body fat monitoring scales is for guidance purposes only. Your physical appearance will not necessarily coincide with the readings for body fat and water, and people of the same age, sex and physical build will not always give the same reading. Consult your doctor to determine the best body water and body fat conditions you should be aiming for.
2. Body fat control scales are not suitable for use by people with the following characteristics:
 - Children under 10 and adults over 80.
 - People with symptoms of fever, oedema or osteoporosis.
 - People undergoing dialysis treatment.
 - People who are taking cardiovascular medication.
 - Pregnant women.
 - Athletes performing over 10 hours of intense exercise per day.
 - People with a heart rate of under 60 beats per minute.

14. INFORMATION FOR THE CORRECT DISPOSAL OF ELECTRICAL AND ELECTRONIC APPLIANCES

At the end of its working life, the product must not be disposed of as urban waste. It must be taken to a special local authority differentiated waste collection centre or to a dealer providing this service. Disposing of a household appliance separately avoids possible negative consequences for the environment

and health deriving from inappropriate disposal and enables the constituent materials to be recovered to obtain significant savings in energy and resources. To remind you that you must collaborate with a selective collection scheme, the symbol shown appears on the product warning you not to dispose of it in traditional refuse containers. For further information, contact your local authority or the shop where you bought the product.

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