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## IRL100 RED LIGHT LAMP Operating Instructions

### Cautions and warnings

- 1.) Read the complete operating instructions before using the IRL 100 Red Light Lamp for the first time.
- 2.) The lamp gets hot in use. Please note the following:
  - Do not touch the lamp while using it or immediately afterwards.
  - Keep the lamp at a safe distance of at least 50cm from any flammable material or object. Do not cover the lamp with textile or other material.
  - Do not look directly in the infrared light emitted by the lamp. When treating the face, keep your eyes closed. Do not sleep.
- 3.) Set up the lamp only on a horizontal surface. The lamp is not made for hanging from a ceiling or wall.
- 4.) Do not expose the lamp to shock or impact. Do not force the lamp to cool down quickly directly after use.
- 5.) Do not use the lamp near water, such as beside a water basin, sink or tub filled with water. Protect the lamp from splashing water.
- 6.) Make sure that children cannot use the lamp. Never leave the lamp unattended when it is connected to power.
- 7.) If in doubt, consult your doctor who will tell you if you can expose your skin to the light of the lamp. This may not be advised if you are taking pain-killing drugs.

### General information

The lamp emits intensive infrared and heat radiation, which penetrates deeply into the skin causing the blood vessels to widen and the circulation to become more vigorous. In this state, the skin and the tissue under the skin are provided more efficiently with the substances they need. This makes the IRL 100 lamp also suitable for cosmetic treatment. Remove all make-up before starting treatment. The radiated heat opens the skin pores so that the skin is cleaned in depth and skin care products, e.g., cream, is more easily taken up by the skin.

### How and how long you should use the lamp

Before starting treatment, make sure that the voltage on the lamp's nameplate is identical with the voltage of the power source you use. The lamp is switched on by plugging it in the socket outlet. Direct the radiation to where you need it with the tilt mechanism of the lamp. Generally, keep a distance of 25 cm to 50 cm between the lamp and the skin you expose to radiation so that you feel the radiation is comfortable. You can vary the duration of radiation, generally the time span is between 15 and 30 minutes.

### Cleaning and maintenance

Allow the lamp to cool slowly after use. Do not force cooling. Clean the lamp after pulling the plug and allowing the lamp to cool. Use a soft, moist (distilled water or alcohol) cloth. When you replace the bulb, make sure that the replacement bulb is of the same rating (100W). The mains lead should only be replaced by an authorized repair center.

### WARRANTY

Your IRL 100 Red Light Lamp was tested for proper function before shipment from the factory. If you return it for repair or inspection, please enclose the bill of sale or cash receipt and a note describing the difficulty you have had with the product or the required testing and mail it to the address below.

In case of defects of manufacture or material, we will repair the product free of charge for a period of 3 years from the date of purchase (date of cash receipt).

The warranty does not apply if the product is used contrary to instructions, handled, stored or cleaned improperly, or as a consequence of any other unusual occurrence. The bulb is a quickly wearing item which is not covered by this warranty. Any further claim, including payment of damages, is excluded.

If you have reason to obtain service under the terms of this warranty, please mail the product with the warranty document to our customer service at the following service addresses:

#### England:

International Business Products Limited  
6 Newfield Court  
West Houghton  
Bolton BL5 3 SA England  
Phone: +44 (0)1942 819695

#### India:

sano health care products pvt. Ltd. **ibp** -Group  
Corporate Office: 28 Landsdown Terrace 1st Floor  
Kolkata 700 026  
Phone: +91 33 2465 0383-85

#### What does this symbol on your product mean?

When it reaches the end of its life cycle, this product must not be disposed of in regular domestic waste but must be taken to a collection point for recycling electrical and electronic appliances (disposal point for used appliances).

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