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**EN** Electronic babyscale ..... p. 04

**DE** Elektronische Babywaage ..... p. 06

**IT** Pesabambini elettronica ..... p. 08

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**NL** Elektronische babyweegschaal ..... p. 12

**Petit**  
Terraillon® 

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Dear customer, Thank for buying this Terraillon *electronic babyscale*.  
We hope you enjoy using it and to ensure your complete satisfaction, we recommend you carefully read this leaflet before use.

## Electronic babyscale

Your Terraillon electronic babyscale has been developed according to quality and safety standards. Its measurement freeze function allows the weight to be displayed even if the baby keeps moving. With a particularly ergonomic design offering comfort and stability, it can memorise a reading and calculate the weight difference between 2 weighs.

### DESCRIPTION

- 1- Memorisation of weight (**M+**)
- 2- Recall of memorised weight (**MRC**)
- 3- Calculates the weight difference between 2 weighs (**▲**)
- 4- Large 78.5 x 37 mm liquid crystal display
- 5- On/off switch (**ON/OFF**)
- 6- Automatic tare (**TARE**)

### FEATURES

Maximum weight: 15 kg  
 Graduation : 10 g from 0 to 6 kg and 15 g from 6 to 15 kg  
 Measurement freeze function  
 Memory Function  
 Indicates the weight difference between two weighs  
 Automatic tare  
 Uses one 9V 6LR61 type battery (not included)  
 Battery indicator

### HOW TO OPERATE YOUR BABYSCALE

#### 1- PUTTING IN THE BATTERY

The appliance uses one 9 V alkaline battery (6LR61 type). Insert the 9V battery into the compartment underneath the babyscale. Respect the polarities +/-.

#### 2- USING THE BABYSCALE

Press the **ON/OFF** button. The screen displays "0g". Place the baby in the scale. The weight is displayed. The measurement freeze function allows a weight to be displayed even if the baby keeps moving.

##### a- Automatic Tare Function

The automatic tare function allows you to weigh baby comfortably lying on a towel without including the weight of the towel.

- Switch on the scales by pressing the **ON/OFF** button.
- Put the towel on the scales.
- Press the **TARE** button. The screen displays "t" then "0g"

- Put your baby comfortably on the towel.

Only the weight of the baby is displayed.

(When you remove baby and towel, a negative weight may appear. This is the corresponding weight of the towel in negative).

##### b- Memorisation of a weight

Once the weight of the baby appears on the screen, you can record it in the memory by pressing the "**M+**" key.

#### c- Recall of recorded weight

By pressing the "MRC" button, the previous weight memorised will appear on the screen with the sign "⊕" to the right. The memorised weight will continue to be displayed for a few seconds.

The recorded weight will remain in the memory until a new weight is recorded or the battery stops working.

#### d- Weight difference between 2 weights

By pressing the "▲" button, the difference between the current weight and the recorded weight will appear on the screen with the sign "▲". If the current weight is greater than the previous weight, the difference will be positive. If it's the other way around, it will be negative.

#### e- Overload

The following message is displayed if the babyscale is overloaded (weight over 15 kg):

"----" then "----". Remove your child immediately or the babyscale will be damaged.

If you wish to change the measuring units from metric to imperial or vice versa, simply activate the conversion switch located on the under-side of the scale.

### RECOMMENDATIONS

- Clean your babyscale with a damp cloth without letting water get into the appliance.
- Remove the 9V battery if the babyscale is not going to be used for a long period.

**Do not throw old batteries away. Place them in a container specially provided for worn batteries so that they may be collected and recycled.**

### GUARANTEE

Your Terraillon electronic babyscale is guaranteed for 2 years against material and production defects. During this period, such defects will be repaired free of charge (proof of purchase must be produced in the case of a claim under guarantee). This guarantee does not cover damages resulting from accidents, incorrect use or negligence. In the case of a claim, please initially contact the shop where you bought your scales.

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**HANSON (UK) LTD**
**TERRAILLON ASIA PACIFIC LTD**